
Notes: To adjust the height of the rod, turn the adjuster assembly counterclockwise, always holding the upper rod.
Slide the stem to the desired height, then turn the adjuster assembly clockwise to lock the upper stem.
Depending on the product model, the pressure gauge can rotate 360° (degrees). For this model, after fitting the manometer to the upper rod, the handle must be screwed into the manometer box/base to fix it.
Max height: 1.55m
Minimum height: 1.00 m


5) Operation / Conditions of Use
The blood pressure check can be affected by the patient's position and respective physiological condition. Before starting the check, make sure that the patient rests for at least 5 minutes, that he has his back and feet supported and that he is not cross-legged. Gently support the patient's forearm and keep the arm at heart level. The procedure needs to be performed in a quiet environment and the patient should not talk or make sudden movements at this time. Failure to observe these recommendations could result in incorrect blood pressure measurements.
Select cuff dimensions suitable for the patient's arm circumference. The applicable dimension is printed on each cuff.
Wrap the cuff around the patient's arm so that the arrow indicating the artery at the bottom of the cuff is over the brachial artery, ie approximately 2.5 cm above the elbow.
The pressure gauge pointer should be indicating "0" – zero.
Place the Stethoscope chestpiece on the artery as indicated on the cuff. Close the pear valve and rapidly inflate to a level about 30 mmHg above the patient's mean systolic pressure.
Carefully open the valve partially to allow deflation at a rate of 2 to 3 mmHg per second.
The first pulsation sound you hear corresponds to systolic (maximum) pressure, continue deflating the cuff by looking closely at the manometer and listening through the Stethoscope for a change in sound intensity until you no longer hear the heartbeat. This is the diastolic (minimum) pressure. Record patient data if necessary.
Open the entire air outlet valve and remove the cuff from the patient.
This product is generally used in beds, hospitals, health centers, domestic use by trained personnel and/or in any circumstance where it is necessary to check blood pressure in human beings.
The product is also used in external transport, such as ambulances, in case of aid and emergencies.
Country legislation must be respected.
6) Operator position in normal use
(Sitting or lying down – comfortable position) forearm up

Auscultation Recommendations
Auscultation of Children 3-12 years: (Stage IV or K4)
Sounds become suddenly and distinctly muffled.
Adult Auscultation: (Phase V or K5)
The sounds cease completely because the artery is no longer compressed and the flow becomes laminar. The pressure indicated on the manometer corresponds to the Diastolic Pressure.
Auscultation of Pregnant Women: (Phase V or K5)
The sounds cease completely because the artery is no longer compressed and the flow becomes laminar. The pressure indicated on the manometer corresponds to the Diastolic Pressure. If sounds are audible with the cuff deflated, stage IV or K4 should be used.

Source: American Heart Association and VI Brazilian Guidelines on Hypertension.

Markings on the Sphygmomanometer

7) Specifications / Regulations
The BIC Sphygmomanometer has an accuracy of ± 3 mmHg at the initial check.
This product maintains the specified safety and performance characteristics at temperatures from 0°C to 40°C with relative humidity not exceeding 85%, non-condensing.
Storage and transport conditions: Store the Sphygmomanometer in a dry and ventilated place, out of reach of children, preferably in its original packaging. Do not bend the tubes.
Transport the product in its original packaging, avoiding falls, knocks, impacts. Respect the conditions stated in the box.
INMETRO/DIMEL Ordinance No. 341/21 - Metrological Technical Regulation
American Standard ANSI/AAMI, SP10: 2002 - Manual Sphygmomanometers
ISO 81060-1: 2012 - Non-invasive sphygmomanometers - Part 1: Requirements and test methods for the non-automated measurement type (ISO 81060-1:2007)
EN ISO 81060-1:2012- Non-invasive sphygmomanometers - Part 1: Requirements and test methods for non-automated measurement type (ISO 81060-1:2007)
Directive 93/42 / EEC amended by Directive 2007/47 / EC
Directive 93/42/EEC amended by Directive 2007/47/EC

8) Reuse: Cleaning, Sterilization and Disinfection
* It is recommended to wash only the cuff with water and neutral soap or cloth moistened with 70° alcohol and dry in the open air. The other parts, cloth slightly dampened with water or liquid alcohol 70°.
*This device cannot be sterilized by the autoclave process.
The use of corrosive products, detergents or polishes is prohibited.
Do not iron.
9) Warranty / Maintenance
Your product, when new, is guaranteed against original manufacturing defects. The warranty period starts from the date of purchase of the product and issuance of the invoice or tax receipt.
CBEMED's obligation is limited to repairing or replacing components considered to be manufacturing defects within the warranty period. These warranties apply to the original purchaser and may not be distributed or transferred to any third party. This warranty does not apply to any damages or failures of the products considered by CBEMED in its analysis to be due to misuse, or use in disagreement with the manufacturer's specifications, accidents including transport damage, falls, malpractice, inadequate maintenance, or outside of the authorized network, tampering, modification, absence of subsequence calibration checks, use of parts or maintenance by unauthorized persons. In addition to voiding the warranty, product repair by unauthorized persons may lead to improper repairs, causing risks to the patient, such as divergent results. Do not allow repairs by unauthorized persons.

Repairs carried out within the warranty period do not extend the warranty period.
These express warranties take precedence over any and all other express or implied warranties, including warranties of merchantability and fitness for a particular purpose, and no other person is authorized to assume on behalf of CBEMED any other warranty in connection with the sale of the product. CBEMED is not responsible for any loss or damage, whether rights, incidental or consequential, that result from the breach of any express warranty, except as set forth herein.
CBEMED covers manufacturing defects free of charge, during the warranty period, as long as the owner presents the purchase receipt or receipt, together with the documents pertaining to the subsequent verifications carried out annually by the RBMLQ-I – Brazilian Network of Legal and Legal Metrology. Quality – Inmetro. This verification is necessary for warranty maintenance, therefore, file such documents as it will be requested.

If the pressure gauge is dropped or has a strong impact, take it to a new check immediately.
CBEMED is responsible for the initial verification carried out at the factory by CBEMED itself and also by the RBMLQ-I – Brazilian Network of Legal Metrology and Quality – Inmetro, in which it approves the product.
The product does not require any preventive maintenance from the operator, other than the annual check of the pressure gauge and clamp, as mentioned above. Consult the RBMLQ-I for your region.
For greater agility and in order to avoid undue freight, when carrying out more careful analyses, CBEMED indicates the shipment of the product by the customer and when the manufacturing defect is confirmed, CBEMED will issue a refund upon proof of value. And with the same objective of speeding up the service, during the case, CBEMED may request evidence of the problem (customer report), such as photos, videos and relevant documents, batch or serial number, always seeking the best and directed service to the customer. client.

10) Warranty Period
Pressure gauge: 5 (years) warranty against manufacturing defects, provided that the circumstances described in this manual are met.
Pear with Valve and Repair: 1 (year) warranty against manufacturing defects.
Clamp: 1 (year) warranty against manufacturing defects.
Cuff: 1 (year) warranty against manufacturing defects.
Rod/Tube: 1 (year) warranty against manufacturing defects.
Star for Pedestal: 1 (year) warranty against manufacturing defects.
Gauge Swivel Part: 1 (year) warranty days against manufacturing defect
If any defect is identified by the purchaser of the product, for the coverage of this warranty term, within the aforementioned terms and conditions, a copy of the purchase invoice or invoice must be sent to CBEMED, otherwise the warranty will not be covered.
Warranty on individual sales/accessories, consult the CBEMED website (www.cbemed.com.br).

11) Environmental Protection
Product not dangerous to the environment. It does not require special care for storage and transport. It has no mercury.
The product contains recyclable materials.
Proper Disposal: Do not dispose of in regular trash. Consult local legislation.

Validity: Undetermined.

MANUFACTURING ACCORDING TO ORDINANCE NO. 46 OF JANUARY 22, 2016.
REGISTRATION WITH THE MINISTRY OF HEALTH No. 80540449001, 80540440007

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Validade: Indeterminada.

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REGISTRO NO MINISTÉRIO DA SAÚDE Nº 80540449001, 80540440007

CONGRATULATIONS ON THE PURCHASE OF YOUR
NEW BIC SPHYGMOMANOMETER

The BIC aneroid sphygmomanometer was subjected to normal and rigorous conditions of use, in addition to tests carried out in laboratories, which prove the quality to guarantee adequate and accurate information in its use. After all, there are decades of tradition in the market of the best sphygmomanometers manufactured in Brazil. Blood pressure measurements determined with this device are equivalent to measurements carried out by a qualified observer using the auscultation method within the limits prescribed by INMETRO Ordinance 341/21.

1) Introduction / Intended Medical Indication

BIC Aneroid Sphygmomanometers are used by healthcare professionals and individuals with immediate auscultation technical training to determine systolic and diastolic blood pressure in humans. Intended Patient Population: Hypertensive or in all cases requiring blood pressure checks in humans.

2) Notices / Usage Parameters

In this manual, a warning identifies a condition or practice that, if not corrected or discontinued immediately, could lead to patient injury, illness, or death.

- Temperature and Environment: Keep in a dry and ventilated place. See item "Specification".

- Pressure: the manometer has a capacity of 0 - 300 mmHg (millimeters of mercury).

- Humidity: Do not store in a wet environment. When cleaning the cuff, allow it to dry completely before storing. Do not wet the gauge.

- Time Limit: Do not leave the cuff on the patient for more than 10 minutes when inflated above 10 mmHg. This can cause disturbances in the patient, disturbing his blood circulation and causing damage to the peripheral nerves.

- Only use BIC brand clamps and accessories. Its replacement by unauthorized (non-original) parts may lead to measurement errors as well as the loss of warranty.

- Avoid falls and excessive heat. Store the device in a cool, dry, ventilated place, out of the reach of children.

- If there is any abnormality or doubts with the equipment, please contact CBEMED immediately. On our website, we also provide our Authorized Technical Assistance. Present this Instruction Manual/Certificate together with the purchase receipt or receipt, in case it is necessary to take your device to a technical assistance.

- Never inflate the cuff above 300 mmHg, you will damage your Sphygmomanometer and VOID THE WARRANTY.

- Keep the Sphygmomanometer, protecting it inside its original packaging, as this also contains the labeling and important data such as batch and product manufacturing.

ATTENTION: It is the responsibility of the instrument holder to carry out the periodic/subsequent verification (calibration) annually with the RBMLQ-1 - Brazilian Network of Legal Metrology and Quality - Inmetro, according to item 5.3.4 of Ordinance No. 341/21.

The product does not require any installation process. Aneroid sphygmomanometers have different dimensions. For each cuff size, it is recommended to use the size compatible with the patient's arm dimension, according to the measurements expressed in the "Components" section.

3) Package Contents / Components

QTY	COMPONENT	TRADITIONAL	ROTATION / TURN	TABLE / WALL
1	Manometer	✓	✓	✓
1	*Cuff/Clamp *(size according to model purchased)	✓	✓	✓
1	Pear with Valve	✓		✓
2	Rods/Tubes		✓	
1	Star for Pedestal		✓	
1	* Knob/screw *(only for the 360° swivel pressure gauge model)		✓	
1	Screw with bushing for Table or Wall model			✓

- Manometer: Metal or plastic box, when using a large manometer, containing aluminum parts. Instrument used to check blood pressure;
- Cuff: Made of hypoallergenic material in cotton or nylon fabric with metal closure or contact closure;
- Cuff: Bag with tubes interconnected between pressure gauges, valve and bulb, made of PVC or vulcanized rubber (depending on the model);
- New Innova Cuff: Made of antiallergic material, free from the inner cuff bag, free from side seams, in nylon fabric and contact closure;
- Pear: Made of high-resistance PVC or vulcanized rubber (depending on the model), used to inflate the cuff;
- Valve: Made of chromed metal, used to control the inlet and outlet of air.
- Rods: Made of treated metal;
- Star for Pedestal: Polypropylene, Nylon and Metal;
- Handle: Polypropylene and Metal.

The models of the BIC brand cuffs are presented in the following sizes: According to Table 1

- Pear: Made of high-resistance PVC or vulcanized rubber (depending on the model), used to inflate the cuff;
- Valve: Made of chromed metal, used to control the inlet and outlet of air.
- Rods: Made of treated metal;
- Star for Pedestal: Polypropylene, Nylon and Metal;
- Handle: Polypropylene and Metal.

The models of the BIC brand cuffs are presented in the following sizes: According to Table 1


*New Innova Adult Cuff: It is 62 cm long, providing greater grip around the circumference, free of internal cuff pouch and side seams.

Scale: 0 to 300 mmHg
Accuracy: ± 3 mmHg at initial check
Resolution: 2 mmHg

4) Caster Assembly

Illustrative images

Step 1: Screw the lower rod into the caster clockwise until the end of the thread. The lower stem contains the regulator assembly.

2nd Step: Insert the upper rod inside the lower rod and adjust the desired height. Then turn the adjuster assembly clockwise to lock.